

HTM2010 Porous Load - Quarterly Tests

| Header | | | |
|-----------------------------|--|-----------------------|----------------------------|
| Authority | Example Health Authority | | |
| Hospital | Example Hospital | | |
| Department | CSSD | | |
| Date of Report | 07/09/2004 | Dates of Tests | 07/09/2004 |
| Sterilizer | Manufacturer Unknown | | |
| | Serial Number No. 01 /014567 1995 | | |
| | Plant reference number None | | |
| | Model Model 2 | | |
| | Usable Chamber Space | N/A litres | Date of Manufacture |
| Service Provider | Isopharm Systems Ltd | | |
| Validation File Ref. | | Data File Ref. | |

| Quarterly Safety Checks | | |
|--|--------------|-------|
| Safety Check | Completed | Notes |
| Check condition of Door Seal: | Yes | |
| Check condition of Closure Mechanism and Limit Switches correctly set: | Yes | |
| Check that the Door is fully engaged by clamps in the locked position: | Yes | |
| Check that the Door Steam Safety interlock is working correctly: | Yes | |
| Check that the Door Safety Edge is working correctly: | Yes | |
| Check that the Door Pressure interlock is working correctly: | Yes | |
| Check that the Chamber Safety Valve is free and discharges at correct press.: | Yes | |
| Check that the Jacket Safety Valve is free and discharges at correct pressure: | Yes | |
| Check condition of Services and test for correct functioning of failure modes: | Satisfactory | |
| Additional checks listed by Competent Person regarding the Pressure Vessel: | Satisfactory | |
| | | |
| | | |

Record Test as PASS or FAIL **PASS**

| Vacuum Leak Test | | | |
|--|---------------------------------------|--------------------------------------|------------------------------------|
| | Before insertion of test Instruments. | After insertion of test Instruments. | After removal of test Instruments. |
| Cycle Counter Number | 328 | 329 | 331 |
| Start time | 00.15.00 hh.mm.ss | 00.15.12 hh.mm.ss | 00.15.15 hh.mm.ss |
| Absolute Gauge Reading (Pump stopped) | 25 mbar | 28 mbar | 30 mbar |
| Absolute Gauge Reading after 5 minutes | 30 mbar | 33 mbar | 35 mbar |
| Absolute Gauge after further 10 minutes | 35 mbar | 38 mbar | 40 mbar |
| Rise in Pressure between Readings | 5 mbar | 5 mbar | 5 mbar |
| Vacuum Leak Rate | 1 mbar/min | 1 mbar/min | 1 mbar/min |
| State PASS/FAIL | PASS | PASS | PASS |

HTM2010 Porous Load - Quarterly Tests

| Automatic Control Test | | | | | | | | | | | |
|--|---------------------------------------|---------------------|--|-------------------------|----------------------------------|---------------------------|-------------------|------------------|---------------|-------------------|------------------|
| | | Cycle No. | 000330 | | | Start time | 12:32:39 hh.mm.ss | | | | |
| Air Removal Stage | Vacuum Switch Point | N/A bar | | | Time to Vacuum Switch Point | | | mins.secs | | | |
| | Pre-Vacuum | Vacuum Hold Time | | | Total Vacuum Duration | | | mins.secs | | | |
| Pulsing | Pulse | | Start | | | Peak | | | Finish | | |
| | | Temp (°C) | Press mB | Time | Temp (°C) | Press mB | Time | Temp (°C) | Press mB | Time | |
| | 1 | 49.2 | -25 | 00:00:06 | 43.1 | -980 | 00:03:11 | 92.7 | -78 | 00:03:32 | |
| | 2 | 92.7 | -78 | 00:03:32 | 73.1 | -756 | 00:03:50 | 94.5 | -74 | 00:04:04 | |
| | 3 | 94.5 | -74 | 00:04:04 | 74.2 | -749 | 00:04:22 | 94.7 | -182 | 00:04:34 | |
| | 4 | 94.7 | -182 | 00:04:34 | 43.1 | -750 | 00:04:52 | 92.6 | -148 | 00:05:05 | |
| | 5 | 92.6 | -148 | 00:05:05 | 74.6 | -735 | 00:05:21 | 92.2 | -142 | 00:05:33 | |
| | 6 | 92.2 | -25 | 00:05:33 | 73.0 | -746 | 00:05:49 | 122.4 | 1452 | 00:06:31 | |
| | 7 | 122.4 | 1452 | 00:06:31 | 112.9 | 407 | 00:06:37 | 120.8 | 1435 | 00:06:52 | |
| | 8 | 120.8 | 1435 | 00:06:52 | 113.6 | 416 | 00:06:58 | 121.3 | 1482 | 00:07:13 | |
| | 9 | 121.3 | 1482 | 00:07:13 | 114.3 | 471 | 00:07:19 | | Not Set | | |
| | 10 | | | | | | | | | | |
| | 11 | | | | | | | | | | |
| | 12 | | | | | | | | | | |
| | 13 | | | | | | | | | | |
| | 14 | | | | | | | | | | |
| | 15 | | | | | | | | | | |
| | | Duration of Pulsing | 0:07:07 | h:mm:ss | Number of Pulses | | 9 | | | | |
| Steam Admission Stage | Time to Attain Sterilizing | 00:31 | mins.secs | | Sterilization temp (ST) selected | | | 134.0 °C | | | |
| Sterilizing Stage | | | | | Start | 1 minute(s) | 2 minutes | 3 minutes | | | |
| | Indicated Pressure (barA) | | | | 3.100 | 3.180 | 3.220 | 3.230 | | | |
| | Recorded Pressure (barA) (Digital**) | | | | 3.060 | 3.190 | 3.220 | 3.220 | | | |
| | Recorded Pressure (barA) (Analogue*) | | | | | | | | | | |
| | Test Pressure Bar | | | | 2.0149 | 2.1844 | 2.2030 | 2.2109 | | | |
| | Indicated Temperature (°C) | | | | 134.0 | 136.0 | 136.0 | 136.0 | | | |
| | Recorded Temperature (°C) (Digital**) | | | | 134.0 | 135.9 | 136.0 | 136.0 | | | |
| | Recorded Temperature (°C) (Analogue*) | | | | | | | | | | |
| | Test Temperature °C | | | | 134.1 | 135.8 | 136.0 | 136.1 | | | |
| | Hold Time | | 00:03:39 | h.mm.ss | | Equilibration Time | 00:00 | mins.secs | | Stage Time | 00:03:39 h.mm.ss |
| Drying Stage | Time to Switch Point | 1.16 | mins.secs | | Vacuum Switch Point | | N/A bar | | | | |
| | Vacuum Hold Period | 0.31 | mins.secs | | Extended Drying Period | | 1.00 mins.secs | | | | |
| | Maximum Vacuum Attained | -0.9781 | bar | | Stage Time | | 3.45 mins.secs | | | | |
| Air Admission Stage | Time to Atmospheric Pressure | N/A mins.secs | | | | | | | | | |
| Cycle Complete Time | 12:49:26 | hh.mm.ss | | Total Cycle Time | 00:16:47 hh.mm.ss | | | | | | |
| Bowie - Dick Test | PASS | | Cycle meets requirements of Automatic Control Test | | | | | Yes | | | |
| ** Digital - taken from Sterilizer digital data log (where fitted). *Analogue - taken from Sterilizer analogue recorder (where fitted) | | | | | | | | | | | |

HTM2010 Porous Load - Quarterly Tests

Thermometric Test for a Small Load

| Cycle No. 000330 | Start time 12:32:39 hh:mm:ss | Start of Plateau Period 00:07:50 hh:mm:ss | Drain 134.1 °C | | | |
|---|--------------------------------------|--|-----------------------|------------|------------------------------|------------------|
| Holding Time Temperature | | 134.0 °C | | | | |
| Channel | Location | Start | 1 minute(s) | 2 minutes | 3 minutes | Time Above |
| Sensor no. 1 | Drain | 134.1 °C | 135.7 °C | 135.9 °C | 136.0 °C | 00:03:39 h:mm:ss |
| Sensor no. 4 | Load | 134.1 °C | 135.8 °C | 136.0 °C | 136.1 °C | 00:03:39 h:mm:ss |
| Indicated | Drain | 134.0 °C | 135.8 °C | 136.0 °C | 136.1 °C | 00:04:02 h:mm:ss |
| Recorded | Drain | 134.0 °C | 135.9 °C | 136.1 °C | 136.1 °C | 00:04:02 h:mm:ss |
| Equilibration Time | | | 00:00 hh:mm:ss | | Drain Fluctuation | 2.0 °C |
| Max. difference Drain - Free Space 1st min. hold time | | | 2.2 °C | | Load Fluctuation | 2.1 °C |
| Max. difference Drain - Free Space remainder of hold time | | | 0.4 °C | | Max. difference Drain - Load | 0.3 °C |
| Holding Time Pressure | | | | | | |
| Channel | Location | Start | 1 minute(s) | 2 minutes | 3 minutes | |
| Sensor no. 13 | Chamber | 2.0149 bar | 2.1844 bar | 2.2030 bar | 2.2109 bar | |
| Indicated | Chamber | 2.050 bar | 2.220 bar | 2.230 bar | 2.250 bar | |
| Recorded | Chamber | 2.060 bar | 2.230 bar | 2.240 bar | 2.250 bar | |
| Post Vacuum (Drying) Stage | | | | | | |
| Stage Time | 3.45 mm:ss | Pressure at end of vacuum hold | -0.9708 bar | | Load sensibly dry | PASS |
| Total Cycle Time | 00:16:47 hh:mm:ss | Cycle meets requirements of Automatic Control Test | | | | PASS |
| Applicable section of HTM 2010 | Applicable HTM 2010 requirements met | | | | Yes | |

HTM 2010 - 13.7

Air Detector - Function Test

| | | |
|---|-------------------|--|
| Cycle No. | Start time | hh:mm:ss |
| Air Detector Setting to Fail Cycle | mbar/°C | Automatic Controller Triggering Point |
| Pressure/Temperature attained by Air Detector | mbar/°C | Leak Rate for Depression dT = Reject* for Small Load |
| Failed Cycle indicated at end of Test | | mbar/min |

* As determined during small load Air Detector Performance Test

HTM 2010 - 11.60

Bowie - Dick Test for Steam Penetration

| | | | |
|---|-------------------------------------|---|-------------------------------------|
| Cycle No. 000330 | Start time 12:32:39 hh:mm:ss | Complete time 12:49:26 hh:mm:ss | Total time 00:16:47 hh:mm:ss |
| Indicated sterilizing temperature | 135 °C | Indicated sterilizing pressure | 2.250 Bar |
| Sterilizing Hold Time | 00:03:39 mins.secs | Plateau time | 00:03:39 mins.secs |
| Maximum Temperature difference between Centre Pack and Drain during Hold Time | 0.3 °C | | |
| Maximum recorded Pressure during Hold Time | 2.2119 Bar | Type of test pack | |
| Uniform Colour change throughout Indicator | PASS | System indicates that Bowie - Dick has been run | Yes |

HTM 2010 - 13.39

HTM2010 Porous Load - Quarterly Tests

Comments

HTM2010 Porous Load - Quarterly Tests

Declaration Section

DECLARATION OF TEST PERSON (STERILIZERS)

1. All test instruments have current calibration certificates.
2. Calibration of the thermometric test instruments has been checked before and after the thermometric tests
3. The periodic checks and tests have been completed and confirm that the sterilizer is safe to use and that commissioning and performance qualification data collected during validation remain valid.

| | | | |
|--------------|------|-----------|------|
| Test Person: | Name | Signature | Date |
|--------------|------|-----------|------|

DECLARATION OF USER

The sterilizer is fit for use. The next periodic tests are due no later than:

| | | | |
|-------|------|-----------|------|
| User: | Name | Signature | Date |
|-------|------|-----------|------|

Calibration Report

| | | |
|---------------------------|------------------|--|
| Job Ref.: | | Thermocouple: Calibration Set 1 |
| | | Date: 07/09/2004 dd/mm/yyyy Time: 15:50:03 |
| Test Person Russell Baker | Sign: VALID ESIG | Date: 09/09/2004 dd/mm/yyyy Time: 11:13:13 |
| Reviewed Russell Baker | Sign: VALID ESIG | Date: 09/09/2004 dd/mm/yyyy Time: 15:38:06 |

| EQUIPMENT | Logger / Recorder | Temp. Ref. Unit |
|--------------|-------------------|--------------------|
| Manufacturer | Agilent | Isotech |
| Model | 34970A | 907 |
| Serial No. | MY41030455 | 34567 |
| Test House | | N.T.P.L |
| Cert. No. | | GRD25082 |
| Renew Date | | 23/06/05 |
| EQUIPMENT | Temp. Ref. Source | Pressure Ref. Unit |
| Manufacturer | Isotech | SI Instruments |
| Model | 907 | TP1-40 |
| Serial No. | 34567 | 5467 |
| Test House | N.T.P.L | RS |
| Cert. No. | GRD25082 | 254156 |
| Renew Date | 23/06/05 | 26/04/05 |

STABILITY SET-UP

0.2 Degrees per minute for 2 minutes. Allowed deviation from Reference 2.0 degrees.
 Reference stability criteria 0.05 degrees for 1 minute. Report after Calibration for 2 minutes every 15 seconds.
Report max deviation allowed 0.5 degrees.

| | | Programmed | Reference |
|-----------|-----------|------------|-----------|
| Setpoints | Low | 100 °C | 100.22 °C |
| | High | 130 °C | 130.25 °C |
| | 1st Check | 121 °C | 121.17 °C |

LOW Point 100.22 °C STABILITY REPORT

Start Time: 15:54:21 Time when stability requirements met: 16:01:22 Elapsed time: 00:07:22
 Reference change: 0.03 °C Maximum sensor change over last minute: 0.06 °C

| | Chan1 | Chan2 | Chan3 |
|-------------|-------|-------|-------|
| Change (°C) | 0.06 | 0.06 | 0.04 |

LOW Point 100.22 °C QUALIFICATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------|----------------|--------|--------|--------|
| 16:01:22 | 100.22 °C | 100.34 | 100.41 | 100.45 |
| | Deviation (°C) | 0.12 | 0.19 | 0.23 |

LOW Point 100.22 °C POST CALIBRATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------|----------------|--------|--------|--------|
| 16:01:38 | 100.20 °C | 100.20 | 100.20 | 100.22 |
| | Deviation (°C) | 0.00 | 0.00 | 0.02 |
| 16:01:53 | 100.20 °C | 100.18 | 100.19 | 100.20 |
| | Deviation (°C) | - 0.02 | - 0.01 | 0.00 |
| 16:02:08 | 100.19 °C | 100.16 | 100.18 | 100.19 |
| | Deviation (°C) | - 0.03 | - 0.01 | 0.00 |
| 16:02:23 | 100.20 °C | 100.21 | 100.24 | 100.26 |
| | Deviation (°C) | 0.01 | 0.04 | 0.06 |
| 16:02:38 | 100.20 °C | 100.24 | 100.22 | 100.24 |
| | Deviation (°C) | 0.04 | 0.02 | 0.04 |
| 16:02:53 | 100.21 °C | 100.22 | 100.23 | 100.23 |
| | Deviation (°C) | 0.01 | 0.02 | 0.02 |
| 16:03:08 | 100.20 °C | 100.16 | 100.19 | 100.19 |
| | Deviation (°C) | - 0.04 | - 0.01 | - 0.01 |
| | MAX DEV. (°C) | 0.04 | 0.04 | 0.06 |

HIGH Point 130.25 °C STABILITY REPORT

Start Time: 16:03:26 Time when stability requirements met: 16:16:22 Elapsed time: 00:12:57

Calibration Report

Job Ref.: _____ Thermocouple: Calibration Set 1
 Date: 07/09/2004 dd/mm/yyyy Time: 15:50:03
 Test Person Russell Baker Sign:VALID ESIG Date: 09/09/2004 dd/mm/yyyy Time: 11:13:13
 Reviewed Russell Baker Sign:VALID ESIG Date: 09/09/2004 dd/mm/yyyy Time: 15:38:06

Reference change: 0.03 °C

Maximum sensor change over last minute: 0.11 °C

| | Chan1 | Chan2 | Chan3 |
|-------------|-------|-------|-------|
| Change (°C) | 0.10 | 0.11 | 0.08 |

HIGH Point 130.25 °C QUALIFICATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 16:16:22 | 130.25 °C | 130.53 | 130.62 | 130.67 |
| Deviation (°C) | | 0.28 | 0.37 | 0.42 |

HIGH Point 130.25 °C POST CALIBRATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 16:16:38 | 130.24 °C | 130.24 | 130.25 | 130.24 |
| Deviation (°C) | | 0.00 | 0.01 | 0.00 |
| 16:16:53 | 130.25 °C | 130.23 | 130.23 | 130.23 |
| Deviation (°C) | | - 0.02 | - 0.02 | - 0.02 |
| 16:17:08 | 130.26 °C | 130.31 | 130.28 | 130.29 |
| Deviation (°C) | | 0.05 | 0.02 | 0.03 |
| 16:17:23 | 130.24 °C | 130.30 | 130.29 | 130.29 |
| Deviation (°C) | | 0.06 | 0.05 | 0.05 |
| 16:17:38 | 130.22 °C | 130.24 | 130.25 | 130.25 |
| Deviation (°C) | | 0.02 | 0.03 | 0.03 |
| 16:17:53 | 130.23 °C | 130.23 | 130.25 | 130.23 |
| Deviation (°C) | | 0.00 | 0.02 | 0.00 |
| 16:18:08 | 130.23 °C | 130.22 | 130.24 | 130.23 |
| Deviation (°C) | | - 0.01 | 0.01 | 0.00 |
| MAX DEV. (°C) | | 0.06 | 0.05 | 0.05 |

CALIBRATION FACTOR AND OFFSET RESULTS

| Channel Number | Status | High Ref. | High Measured | Low Ref. | Low Measured |
|----------------|----------|-----------|---------------|----------|--------------|
| 1 | IN SPEC. | 130.25 | 130.53 | 100.22 | 100.34 |
| 2 | IN SPEC. | 130.25 | 130.62 | 100.22 | 100.41 |
| 3 | IN SPEC. | 130.25 | 130.67 | 100.22 | 100.45 |

FIRST CALIBRATION CHECK Point 121.17 °C STABILITY REPORT AFTER ADJUSTMENT

Start Time: 16:18:26 Time when stability requirements met: 16:25:12 Elapsed time: 00:06:47
 Reference change: 0.02 °C Maximum sensor change over last minute: 0.10 °C

| | Chan1 | Chan2 | Chan3 |
|-------------|-------|-------|-------|
| Change (°C) | 0.09 | 0.08 | 0.10 |

FIRST CALIBRATION CHECK Point 121.17 °C QUALIFICATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 16:25:12 | 121.17 °C | 121.28 | 121.28 | 121.28 |
| Deviation (°C) | | 0.11 | 0.11 | 0.11 |

FIRST CALIBRATION CHECK Point 121.17 °C REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 16:25:28 | 121.18 °C | 121.22 | 121.24 | 121.24 |
| Deviation (°C) | | 0.04 | 0.06 | 0.06 |
| 16:25:43 | 121.16 °C | 121.21 | 121.22 | 121.22 |
| Deviation (°C) | | 0.05 | 0.06 | 0.06 |
| 16:25:58 | 121.15 °C | 121.22 | 121.23 | 121.23 |
| Deviation (°C) | | 0.07 | 0.08 | 0.08 |
| 16:26:13 | 121.16 °C | 121.27 | 121.27 | 121.26 |
| Deviation (°C) | | 0.11 | 0.11 | 0.10 |

Calibration Report

| | | | | | |
|---------------------------|--|--|--|-----------------|--|
| Job Ref.: | | | | | Thermocouple: Calibration Set 1 |
| | | | | | Date: 07/09/2004 dd/mm/yyyy Time: 15:50:03 |
| Test Person Russell Baker | | | | Sign:VALID ESIG | Date: 09/09/2004 dd/mm/yyyy Time: 11:13:13 |
| Reviewed Russell Baker | | | | Sign:VALID ESIG | Date: 09/09/2004 dd/mm/yyyy Time: 15:38:06 |

| | | | | |
|----------------|-----------|--------|--------|--------|
| 16:26:28 | 121.17 °C | 121.30 | 121.29 | 121.28 |
| Deviation (°C) | | 0.13 | 0.12 | 0.11 |
| 16:26:43 | 121.16 °C | 121.28 | 121.28 | 121.27 |
| Deviation (°C) | | 0.12 | 0.12 | 0.11 |
| 16:26:58 | 121.16 °C | 121.23 | 121.24 | 121.23 |
| Deviation (°C) | | 0.07 | 0.08 | 0.07 |
| MAX DEV. (°C) | | 0.13 | 0.12 | 0.11 |

| | | |
|---------------------------------------|-----------------------------|---|
| Cycle No. 000330 | Isopharm Systems Ltd | 15/09/2004 |
| Serial No. No. 01 /014567 1995 | Example Hospital | |
| Porous Load Demo | | C:\logsys\AC000001\data\000330.mbf |

| | | | |
|---------------------|-------------|-------------|-------------------|
| Test Person: | Sign | Date | <i>dd/mm/yyyy</i> |
| Reviewed By: | Sign | Date | <i>dd/mm/yyyy</i> |

| Test Details | | | |
|-----------------------------|---------------------------|-------------------------------|--|
| Cycle Number 000330 | Comments | | |
| Date 15/09/2004 | | | |
| Start Time: 12:32:39 | End Time: 12:49:26 | Operator Russell Baker | |
| Batch | | Product | |
| Operating Cycle Ref | | User | |

| | |
|-----------------------------|------------------------------|
| Logger | Serial No. MY41030456 |
| Manufacturer Agilent | Model 34970A |
| Cert. No. | Renew Date |
| Test House | |

| Test Specification | General | Porous Load Demo |
|---------------------------------|----------|-------------------------------|
| Sterilising temperature | 134.0 °C | |
| Scan Interval (secs) | 1 | Version: |
| Auto stop logging after: | dd:hh:mm | Date of last revision: |
| Stage Analysis | ON | Log off temperature °C |
| | | Log off channel |

| Test Specification | Lethality | Porous Load Demo |
|----------------------------|-----------|---------------------|
| Base temperature: | 0.0 °C | Z Value: 5.0 |
| Selected Channels | | |
| Start lethality on: | 0.00 °C | |
| End lethality on: | | |


| CHANNEL CONFIGURATION Simulation | | | | | | |
|----------------------------------|--------|---------------------|--------|---------------|------------|--|
| No. | Label: | Channel Type: | Units: | Chart Colour: | Location: | |
| 1 | Tmp 1 | Type T range: 100mV | °C | Light Red | Drain | |
| 4 | Tmp 4 | Type T range: 100mV | °C | Light Blue | Pack | |
| 7 | Tmp 7 | Type T range: 100mV | °C | Light Cyan | Top | |
| 11 | Tmp11 | Type T range: 100mV | °C | Light Magenta | Bottom | |
| 12 | Tmp12 | Type T range: 100mV | °C | Light Green | Free Space | |
| 13 | Press | +/- 10V range: 10V | Bar | Black | Pressure | |

| Calibration | | | | | | | | Calibration Set 1 | |
|-------------|------------|----------|------------|-------------|------------|-------------|-----|-------------------|--|
| No. | Date | Time | Measured 1 | Reference 1 | Measured 2 | Reference 2 | | | |
| 1 | 07/09/2004 | 16:25:14 | 130.53 | 130.25 | 100.34 | 100.22 | °C | | |
| 4 | | | | | | | °C | | |
| 7 | | | | | | | °C | | |
| 11 | | | | | | | °C | | |
| 12 | | | | | | | °C | | |
| 13 | | | | | | | Bar | | |

Cycle No. 000330
 Serial No. No.01 /014567 1995

Isopharm Systems Ltd
 Example Hospital


C:\logsys\AC000001\data\000330.mbf
 15/09/2004
 Porous Load Demo



Isopharm Systems Ltd
 Example Hospital

Cycle No. 000330
 Serial No. No. 01 /014567 1995

15/09/2004
 Porous Load Demo



Calibration Check Report

Job Ref.: _____ Thermocouple: Calibration Set 1
 Date: 15/09/2004 dd/mm/yyyy Time: 12:50:24
 Operator: _____ Sign: _____ Date: _____
 Reviewed By: _____ Sign: _____ Date: _____

| EQUIPMENT | Logger / Recorder | Temp. Ref. Unit |
|--------------|-------------------|--------------------|
| Manufacturer | Agilent | Isotech |
| Model | 34970A | 907 |
| Serial No. | MY41030456 | 34567 |
| Test House | | N.T.P.L |
| Cert. No. | | GRD25082 |
| Renew Date | | 23/06/05 |
| EQUIPMENT | Temp. Ref. Source | Pressure Ref. Unit |
| Manufacturer | Isotech | SI Instruments |
| Model | 907 | TP1-40 |
| Serial No. | 34567 | 5467 |
| Test House | N.T.P.L | RS |
| Cert. No. | GRD25082 | 254156 |
| Renew Date | 23/06/05 | 26/04/05 |

STABILITY SET-UP

0.2 Degrees per minute for 2 minutes. Allowed deviation from Reference 2.0 degrees.
 Reference stability criteria 0.05 degrees for 1 minute. Report after Calibration for 2 minutes every 15 seconds.
 Report max deviation allowed 0.5 degrees.

| Setpoints | 1st Check | Programmed | Reference |
|-----------|-----------|------------|-----------|
| | | 134 °C | 134.02 °C |

FIRST CALIBRATION CHECK Point 134.02 °C STABILITY REPORT AFTER ADJUSTMENT

Start Time: 12:50:31 Time when stability requirements met: 12:52:33 Elapsed time: 00:02:03
 Reference change: 0.02 °C Maximum sensor change over last minute: 0.03 °C

| Change (°C) | Chan1 | Chan2 | Chan3 |
|-------------|-------|-------|-------|
| | 0.03 | 0.03 | 0.03 |

FIRST CALIBRATION CHECK Point 134.02 °C QUALIFICATION REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 12:52:33 | 134.02 °C | 134.06 | 134.02 | 134.00 |
| Deviation (°C) | | 0.04 | 0.00 | - 0.02 |

FIRST CALIBRATION CHECK Point 134.02 °C REPORT

| Time | Ref. | Chan1 | Chan2 | Chan3 |
|----------------|-----------|--------|--------|--------|
| 12:52:49 | 134.02 °C | 134.06 | 134.04 | 133.99 |
| Deviation (°C) | | 0.04 | 0.02 | - 0.03 |
| 12:53:04 | 134.02 °C | 134.08 | 134.03 | 133.98 |
| Deviation (°C) | | 0.06 | 0.01 | - 0.04 |
| 12:53:19 | 134.02 °C | 134.07 | 134.03 | 133.99 |
| Deviation (°C) | | 0.05 | 0.01 | - 0.03 |
| 12:53:34 | 134.02 °C | 134.07 | 134.03 | 133.99 |
| Deviation (°C) | | 0.05 | 0.01 | - 0.03 |
| 12:53:49 | 134.02 °C | 134.06 | 134.04 | 134.00 |
| Deviation (°C) | | 0.04 | 0.02 | - 0.02 |
| 12:54:04 | 134.02 °C | 134.07 | 134.04 | 133.98 |
| Deviation (°C) | | 0.05 | 0.02 | - 0.04 |
| 12:54:19 | 134.02 °C | 134.06 | 134.03 | 134.00 |
| Deviation (°C) | | 0.04 | 0.01 | - 0.02 |
| MAX DEV. (°C) | | 0.06 | 0.02 | - 0.04 |